

## **Ciliary Body (iridocyclectomy specimen)**

### **Special arrangements**

No special arrangements necessary.

### **Fixation**

10% neutral buffered formalin.

### **Special instructions**

- ❖ Because this specimen is usually taken for tumor, a study of the resection margins is important. Please make a sketch on the Ocular Pathology Consultation Request form of the location of the tumor in the eye so that the medial and temporal resection margins can be distinguished in the specimen.
- ❖ Follow the instructions listed under General Information for Routine Specimen Collection for proper packaging and labeling of routine specimens and how to fill out the Ocular Pathology request form.
- ❖ Follow the procedure under General Information on Collection and Handling of Specimens for proper specimen sending information.